

SEDEX Low Temperature Evaporative Light Scattering Detector

SAMPLE: MONO, DI, TRISACCHARIDES


HPLC
Gradient

Peaks

dR: Deoxyribose
Rh: Rhamnose
F: Fructose
G: Glucose
S: Sucrose
R: Raffinose

Gradient:

Time	% A	% B	% C
0 min	80.0	19.8	0.2
5 min	80.0	19.8	0.2
7 min	45.0	54.8	0.2



CONDITIONS

Column	Zorbax SIL 7 μ 100 x 4 mm
Injected Volume	20 μ l
Eluent	A: CH ₂ Cl ₂ B: Methanol C: Water
Flow Rate	1.5 ml/min
Pressure	2 bars
Gas	Air
Temperature	35° C

Analyzed by: LCBA 45000 Orleans FRANCE

S.E.D.E.R.E.
PARC VOLTA - BP 27
9, rue Parmentier
94141 Alfortville Cedex FRANCE

Application #74