


# SEDEX Low Temperature Evaporative Light Scattering Detector

**SAMPLE: LINOLEIC ACID**

**HPLC**



## CONDITIONS

Column	<i>Zorbax ODS</i> <i>250 x 4.6 mm</i>
Injected Volume	<i>10 µl</i>
Eluent	<i>Methanol/5 % Trifluoroacetic Acid</i> <i>(75/25)</i>
Flow Rate	<i>1.5 ml/min</i>
Pressure	<i>2.0 bars</i>
Gas	<i>Helium</i>
Temperature	<i>35° C</i>

Analyzed by: CE Saclay Molecules Marquees 91191 Gif Sur Yvette FRANCE