

SEDEX Low Temperature Evaporative Light Scattering Detector

SAMPLE: POLYOXYETHYLENE ALCOHOLS

HPLC



The lower temperature used in the SEDEX detector provides a significant advantage as it can detect compounds at a lower temperature. This application shows that many compounds which are readily observed at 35°C are not observable at 65°C (due to thermal instability, thus the analysis can only be performed via the SEDEX detector).

Analysis A Temperature = 30°C
Analysis B Temperature = 65°C

Analysis A can only be performed via the SEDEX Detector

Peaks

R ₁₂ OH	10.5 mg/ml
R ₁₂ (OCH ₂ CH ₂)OH	10.2 mg/ml
R ₁₂ (OCH ₂ CH ₂) ₂ OH	11.8 mg/ml
R ₁₂ (OCH ₂ CH ₂) ₃ OH	9.2 mg/ml
R ₁₂ (OCH ₂ CH ₂) ₄ OH	10.2 mg/ml
R ₁₂ (OCH ₂ CH ₂) ₆ OH	10.4 mg/ml
R ₁₂ (OCH ₂ CH ₂) ₈ OH	9.5 mg/ml
R ₁₂ (OCH ₂ CH ₂) ₉ OH	13.0 mg/ml



CONDITIONS

Column	Lichrospher NH ₂ 250 x 4.6 mm
Eluent	Hexane/iso-propyl alcohol/water
Flow Rate	1.0 ml/min
Pressure	2 bars
Gas	Nitrogen

S.E.D.E.R.E.

PARC VOLTA - BP 27

9, rue Parmentier

94141 Alfortville Cedex FRANCE

Application #89