


SEDEX Low Temperature Evaporative Light Scattering Detector

SAMPLE: UREA

Micro HPLC

Gradient

- 1) Dimyristoylphosphatidyl glycerol, 1.2 µg/µl
- 2) Mannitol, 14 µg/µl
- 3) Dimyristoylphosphatidyl choline, 2.8 µg/µl



CONDITIONS

Column	<i>Nucleosil Diol, 7 µm 250 x 4.0 mm</i>
Injected sample	<i>20 µl</i>
Eluent	<i>i-Propanol/Methanol/H₂O/NH₄OH (con) (57.5/36.5/5.9/0.1)</i>
Flow Rate	<i>1.2 ml/min</i>
Pressure	<i>2.2 bars</i>
Gas	<i>Nitrogen</i>
Temperature	<i>40° C</i>

Analyzed By : Richard Scientific, Novato CA 94949 USA

S.E.D.E.R.E.

PARC VOLTA - BP 27

9, rue Parmentier

94141 Alfortville Cedex FRANCE

Application #95