


SEDEX Low Temperature Evaporative Light Scattering Detector

SAMPLE: FAT SOLUBLE VITAMIN STANDARDS

HPLC

- 1 = Vitamin A
- 2 = Vitamine E
- 3 = Vitamin K
- 4 = β -Carotene



CONDITIONS

Column	<i>Kromasil ODS, 5 μ (250 x 4.6 mm)</i>
Injected sample	<i>20 μl (~ 5 ppm Vitamins E,K, β-carotene, Vitamin D unknown)</i>
Eluent	<i>Acetonitrile/THF (80/20)</i>
Flow Rate	<i>1 ml/min</i>
Pressure	<i>2.2 bars</i>
Gas	<i>Nitrogen</i>
Temperature	<i>40° C</i>

Analyzed By : Richard Scientific, Novato CA 94949 USA

S.E.D.E.R.E.
PARC VOLTA - BP 27
9, rue Parmentier
94141 Alfortville Cedex FRANCE

Application #100