


# SEDEX Low Temperature Evaporative Light Scattering Detector

**SAMPLE : HYDROLIZED SOYBEAN OIL**

**HPLC**



1. Triglycerides
2. Diglycerides
3. Monoglycerides

## CONDITIONS

Column	<i>Nucleosil 70H (diol), 250 x 4 mm</i>
Injected volume	<i>10 µl 1 mg/ml</i>
mobile phase	<i>85 % isooctane ; 10 % IPA ; 5 % MeOH</i>
Flow rate	<i>1 ml/min</i>
Pressure	<i>2.3 bars</i>
Gas	<i>Nitrogen</i>
Gain	<i>8</i>
Temperature	<i>40°C</i>

Analyzed by : Richard SCIENTIFIC U.S.A.