


# SEDEX Low Temperature Evaporative Light Scattering Detector

SAMPLE: ESTRADIOL AND ESTRONE

HPLC



## CONDITIONS

|                 |   |
|-----------------|---|
| Column          | <i>Zorbax Sil - 5 <math>\mu</math>m</i><br><i>25 x 4.6 mm</i> |
| Injected Volume | <i>20 <math>\mu</math>l</i>                                   |
| Eluent          | <i>Cyclohexane/Ethyl Acetate (95/5)</i>                       |
| Gas             | <i>Nitrogen</i>   |
| Temperature     | <i>35° C</i>  |